




Maternity COVID TESTING QUICK TIPS

	Roche- PCR	Abbott- Rapid	Community
When to use	<ul style="list-style-type: none"> Admitted (not triage) unscheduled patients presenting in early labor Same day inductions Admitted antepartum patients 	<ul style="list-style-type: none"> Unscheduled patients presenting in active labor 	<ul style="list-style-type: none"> Elective C-Section patients 72 hours from admission Elective Inductions 72 hours from admission Patients scheduled within 24-48 hours of admission
Order number	SARS-CoV-2 PCR LAB98416	SARS-CoV-2 Rapid Test LAB98414	SARS-COV-2 LAB98415
Collection type	Nasopharyngeal	Nasal	Nasopharyngeal
Swab type	Dry Swab- sterile Miraclean swab 	Foam Puritan swab 	Floq nasopharyngeal swab with medium 
How to send	Follow campus specific process *** do not tube	Follow campus specific process *** do not tube	Follow campus specific process *** do not tube FSED- Courier to Lab
Turn-around time	Up to 24 hours	15-20 minutes	24-48 hours- Positive Results go to the NOCC and then called to the Physician